

# Gentech<sup>®</sup> Generators

## Molybdenum - 99Mo Technetium - 99mTc - Sterile Generator

[Product Information \(PI\)](#)

[Consumer Medical Information \(CMI\)](#)

[Directions for Use](#)

[Generator Returns – Australian Customers](#)

[Generator Returns – International Customers](#)

[Elution Profile](#)

[Gentech Accessories.pdf](#)

---

Product Code	10000
Pack Size	20GBq, 30GBq, 40GBq, 50GBq, 60GBq, 80GBq, 100GBq, 120GBq, 370GBq
Indications	Sodium pertechnetate [ <sup>99m</sup> Tc] is used for scintigraphy, principally of the brain and thyroid. It can also be used to prepare various technetium-99m labelled injections for selective organ imaging especially of the liver, lung, bone and kidney.
Calibration	0900 hours Sydney time Monday and Thursday
Availability	7 days' notice required
Expiry	14 days from calibration
Aust R	72820 75859